

Gebauer's Single Use PainEase One™

Single-Use Topical Anesthetic Designed for Modern Clinical Practice

Backed by over 20 years of clinical research

Designed for One. Trusted by Many.

Healthcare organizations face growing pressure to improve efficiency, standardize processes, and optimize resources while maintaining high-quality patient care. PainEase One™ is a ready-to-use, single-use topical anesthetic designed to help clinicians and healthcare leaders meet these challenges. For more than a century, healthcare providers have trusted Gebauer Company for innovative solutions that enhance patient comfort.

Why PainEase One™?

Ready-to-Use Convenience

Enables clinical teams to focus on patient care rather than product management

Single-Use Format

Individually packaged for one-time use, supporting standardized protocols and consistent administration practices

Workflow Efficiency

Designed to integrate seamlessly into procedural and outpatient settings to help streamline care delivery



Benefits for Healthcare Teams

- Begins working in 4-10 seconds
- Non-drug, non-flammable
- Ready when needed, easy to use, single-use
- Can be used by any licensed healthcare practitioner without the order of a physician

A Smarter Approach to Topical Anesthetics

PainEase One™ delivers a practical, single-use solution designed to support clinicians, optimize workflows, and help healthcare organizations achieve their clinical and operational objectives.



 **Gebauer Company**
www.gebauer.com

Gebauer's Single Use PainEase One™

Single-Use Topical Anesthetic Designed for Modern Clinical Practice

Backed by over 20 years of clinical research

Indications and Usage

PainEase One™ delivers a practical, single-use solution to temporarily control pain associated with:

- Immunizations
- IV starts
- Joint injections
- Blood draws
- Incision and drainage of small abscesses
- Foreign body removal
- Lancing boils
- Skin tag removal
- Suturing and suture removal
- Cosmetic procedures

Important Risk and Safety Information for Gebauer's PainEase One™

Consult your pediatrician when using on children 4 years old and younger. Do not use on large areas of damaged skin, puncture wounds, animal bites or serious wounds. Do not spray in eyes. Over spraying may cause frostbite. Freezing may alter skin pigmentation. Use caution when using product on persons with poor circulation. Apply only to intact oral mucous membranes. Do not use on genital mucous membranes. The thawing process may be painful and freezing may lower resistance to infection and delay healing. If skin irritation develops, discontinue use. CAUTION: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

PainEase One 0.4 fl. oz. in Stream **P/N 0386-0008-05**

3 Step Usage for PainEase One™



Prep

- Have all the necessary equipment ready
- Prepare the skin per your protocol
- Prepare the patient for the cold effect



Spray

- Hold the can 3 to 7 inches from treatment site, about a can's length away
- Spray steadily 4 to 10 seconds or until the skin begins turning white, whichever comes first
- Do not spray longer than 10 seconds



Perform

- Perform the procedure
- The anesthetic effect lasts about one minute

Refer to application instructions for full details

